




TOGETHER Enteral Safety Rises to a New Level of Protection




BY COMBINING THE NEW MEDFUSION® 3500 V5 SYRINGE INFUSION PUMP AND OUR NEOMED LINE OF ORAL/ENTERAL SYRINGES, A NEW STANDARD OF ENTERAL SAFETY IS ACHIEVED.

The Medfusion® 3500 V5 Pump uses smart pump technology that is designed to detect NeoMed's oral/enteral syringes when placed in the pump's syringe chamber. The Medfusion 3500 V5 Pump can be programmed to contain separate enteral feeding programs (enteral delivery of volume over time). This technology helps prevent errors associated with the misuse of IV syringes when delivering enteral nutrition.

To begin a feeding, the display of the Medfusion 3500 V5 Pump shows  on the screen for verification. Selecting syringe size, volume and time completes the process and feeding begins. Expensive and cumbersome dedicated enteral pumps are no longer required.

For additional information, please contact your local NeoMed distributor or reach NeoMed at www.neomedinc.com




For safe, accurate enteral nutrition delivery that meets JCAHO requirements for enteral feedings

Images provided courtesy of Smiths Medical ASD, Inc., St. Paul, Minnesota

Hospital Product List

Product Code	Product Description	Quantity
Individually Packaged Dispensers (Pump Compatible)		
NM-S1EO**	Individually Packaged Oral / Enteral Dispensers 1ml	200/Box
NM-S3EO*	Individually Packaged Oral / Enteral Dispensers 3ml	200/Box
NM-S6EO	Individually Packaged Oral / Enteral Dispensers 6ml	200/Box
NM-S12EO	Individually Packaged Oral / Enteral Dispensers 12ml	200/Box
NM-S20EO	Individually Packaged Oral / Enteral Dispensers 20ml	100/Box
NM-S35EO	Individually Packaged Oral / Enteral Dispensers 35ml	100/Box
NM-S60EO	Individually Packaged Oral / Enteral Dispensers 60ml	100/Box
NM-S100EO*	Individually Packaged Oral / Enteral Dispensers 100ml	100/Box

Sterile Individually Packaged Syringes (Pump Compatible)

SD-S1EO**	Sterile Oral / Enteral Syringes 1ml	200/Box
SD-S3EO*	Sterile Oral / Enteral Syringes 3ml	200/Box
SD-S6EO	Sterile Oral / Enteral Syringes 6ml	200/Box
SD-S12EO	Sterile Oral / Enteral Syringes 12ml	200/Box
SD-S20EO	Sterile Oral / Enteral Syringes 20ml	100/Box
SD-S35EO	Sterile Oral / Enteral Syringes 35ml	100/Box
SD-S60EO	Sterile Oral / Enteral Syringes 60ml	100/Box
SD-S100EO*	Sterile Oral / Enteral Syringes 100ml	100/Box


* Validation is pending for the NM-S3EO, SD-S3EO, NM-S100EO, and SD-S100EO

** The NM-S1EO and SD-S1EO are not characterized with the Medfusion 3500 Version 5 Pump

NeoMed offers the only complete enteral delivery system that complies with all recommendations set by The Joint Commission Journal on Quality and Patient Safety, May 2008, Vol. 34 Num. 5 pg. 285-292, "Enteral Feeding Misconnections: A Consortium Position Statement"

"Do not modify or adapt IV feeding devices because doing so may compromise the safety features incorporated into their design."

(The Joint Commission Journal on Quality and Patient Safety, Enteral Feeding Misconnections: A Consortium Position Statement, pg. 289)



Orange graduation markings and text on the barrel of the syringe are a registered trademark of NeoMed, Inc.

8-0111