



Case Cart

This unit is constructed of Oak, Cherry or Maple and is available in your choice of custom stain colors and hardware. The Case Cart is available in a variety of door and drawer configurations. Standard features include , 14" pop-up work surface,(2nd pop-up optional), wood or laminated top, Hafele Aximate 270° hinged doors, MedCastor 5" two position locking castors, mounted to an aluminum substrate. Beveled or flat slab solid hardwood drawer fronts. Solid hardwood doors. Dovetail constructed hardwood drawer boxes. Accuride 3832 self-closing, full extension drawer slides. Vinyl T-mold edged base deck.



Dimensions:

32.5"W x 24" D x 34"H

Homewood Health Care Products Inc.
6100 Avery Road
Dublin, Ohio 43016

 Contract Holder

www.homewoodhealthcare.com



(P) 614-766-4862 (P) 800-422-6264 (F) 614-766-7248