



## Fetal Monitor Cart

This unit is constructed of Oak, Cherry or Maple and is available in your choice of custom stain colors and hardware. Designed to fit your Fetal Monitoring System. The Fetal Monitor Cart is available in a variety of door and drawer configurations. Standard features include , flip down CPU drawer front, self closing concealed door hinges. Shock metal keyboard tray pull out, wood or laminated top. Hafele Aximate 270° hinged doors, MedCastor 5" two position locking castors, mounted to an aluminum substrate. Beveled or flat slab solid hardwood drawer fronts. Solid hardwood doors. Dovetail constructed hardwood drawer boxes. Accuride 3832 self-closing, full extension drawer slides. Vinyl T-mold edged base deck. Optional 14" pop-up work surface.



### Dimensions:

24.5"W x 25.5"D x 42"H

Homewood Health Care Products Inc.  
6100 Avery Road  
Dublin, Ohio 43016



[www.homewoodhealthcare.com](http://www.homewoodhealthcare.com)



(P) 614-766-4862 (P) 800-422-6264 (F) 614-766-7248