



Infant Care Center®

This unit is constructed of Oak, Cherry or Maple and is available in your choice of custom stain colors and hardware. The Infant Care Center® is available in a variety of door, drawer and internal shelf configurations. Standard features include; Two doors, Hafele XL-slide pocket door system. Grass 120° concealed hinges. Beveled or flat slab solid hardwood drawer fronts. Solid hardwood doors. Dovetail constructed hardwood drawer boxes. Accuride 3832 self-closing, full extension drawer slides. Includes pull out drop side bassinet with mattress. Choice of Wall Mount Infant Warmer. Medical Gas and Electrical Outlets optional. Neo-Puff Resuscitator, optional.



Dimensions:

Customizable to the space available.

Homewood Health Care Products Inc.
6100 Avery Road
Dublin, Ohio 43016

 Contract Holder

www.homewoodhealthcare.com



(P) 614-766-4862 (P) 800-422-6264 (F) 614-766-7248