



Infant Warmer Cabinet



This unit is constructed of Oak, Cherry or Maple and is available in your choice of custom stain colors and hardware. The Warmer Cabinet is available in a variety of solid hardwood door configurations. Standard features include; Hafele Aximate 270° Bi-Fold doors, standard crown Molding and ventilation grommets on top. Choice of Wall Mount Infant Warmer. Medical Gas and Electrical Outlets optional. Neo-Puff Resuscitator, optional. Dropside or Tub Style Bassinet available from Homewood.

Dimensions:

44" W X 32" H X 14 ½ D

Dimensions vary based on warmer used

Homewood Health Care Products Inc.
6100 Avery Road
Dublin, Ohio 43016

 Contract Holder

www.homewoodhealthcare.com



(P) 614-766-4862 (P) 800-422-6264 (F) 614-766-7248