



## Bassinet, 29" 31" 32" 34"

This unit is constructed of Oak, Cherry or Maple and is available in your choice of custom stain colors and hardware. The Bassinet is available in a variety of door and drawer configurations. Standard features include a pull-out work surface, Lexan chart holder, One 2 ½ " T-block, MedCastor 5" two position locking castors, mounted to an aluminum substrate. Beveled or flat slab solid hardwood drawer fronts. Solid hardwood doors. Dovetail constructed hardwood drawer boxes. Accuride 3832 self-closing, full extension drawer slides. Lexan chart holder and a vinyl T-mold edged base deck. Tub and mattress not included



### Dimensions:

32.5"W x 21"D x 29, 31, 32, 34" H  
(Height to top of mattress).

Homewood Health Care Products Inc.  
6100 Avery Road  
Dublin, Ohio 43016

[www.homewoodhealthcare.com](http://www.homewoodhealthcare.com)

(P) 614-766-4862 (P) 800-422-6264 (F) 614-766-7248