

Kirlin Recessed Motorized Exam Lights

Kirlin Lighting, Innovation Design and Performance Since 1895

www.kirlinlighting.com

Kirlin meets these needs with 5 different styles of Motorized Medical Procedural Lights

- Used in L&D, NICU, PICU, ICU, ER, Cath Labs, Special procedural areas and more.
- Shadow free, 60 degree high angle illumination for optimal visualization of deep body cavities.
- 25,600 foot candles at the site.
- Patented instantaneous and precise control of light direction and intensity with infrared remote and wall control.
- Dichroic (cold mirrored) reflectors transmit 70% less heat, assuring patient and practitioner comfort.
- Inconspicuous and fully recessed, low voltage designs
- Shadow free illumination.
- Complete 3 year warranty.

Easy to install and service. All service is through the light aperture and all components are modular for fast and simple service. Relamping is quick and easy