

Neonatal Resuscitation Masks





Neonatal Resuscitation Masks

- Easy to use
- Easy to seal
- Wide size range
- Extra small sizes for premature neonates

Product Code	Sizes in diameter
RD803-10	35mm (1.38in)
RD804-10	42mm (1.65in)
RD805-10	50mm (1.97in)
RD806-10	60mm (2.36in)
RD807-10	72mm (2.83in)

The product does not contain Latex or DEHP



INTERNATIONAL P O Box 14 348, Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web Site: www.fphcare.com
USA Tel: 1800 446 3908 or +1 949 470 3900 Fax: +1 949 470 3933 **AUSTRALIA** Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **AUSTRIA** Tel: 0800 29 31 23 Fax: 0800 29 31 22
BENELUX Tel: +31 78 644 0924 Fax: +31 78 644 0914 **CHINA** Tel: +86 20 3205 3486 Fax: +86 20 3205 3495 **FRANCE** Tel: +33 1 6446 5201 Fax: +33 1 6446 5221
GERMANY Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **INDIA** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **IRISH REPUBLIC** Tel: 1800 409 011 **ITALY** Tel: +39 06 7839 2939
 Fax: +39 06 7814 7709 **JAPAN** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **MIDDLE EAST** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **SPAIN** Tel: +34 902 013 346 Fax: +34 902 013 379
SWITZERLAND Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)** Tel: +44 1628 626 136 Fax: +44 1628 626 146 **Nthn Ireland** Tel: 0800 132 189